

DF-330

SURGICAL & GYNECOLOGY SUCTION UNIT (Aspirator)



www.doctorsfriend.com.tw



DOCTOR'S FRIEND

CE 2460 ISO 13485

DF-330

SURGICAL & GYNECOLOGY SUCTION UNIT (Aspirator)



The SURGICAL & GYNECOLOGY SUCTION UNIT has been designed for professional aspiration of body fluids, tissues or bones of the patient during or after surgery with the latest innovation in high-vacuum suction. Two-in-one function machine improve work efficiency.

SURGICAL & OBS. suction function may be switched by selective knob; **small** and **large calibers** are equipped in two collection bottles for general surgery or OBS.



A spacious storage platform for instrument.



The tube holder is to hold the suction tube, catheter or suction instrument.



The aluminum bottle rack offers more space for collection bottles and enlarges the capacity. According to user's preference, reusable and disposable bottles are able to install.

Autoclavable and reusable PC collection bottle is factory equipped.



The fan is designed to prevent heat and make pump works efficiency.

Lower center of gravity makes it moving stable and smoothly; caster with locking brake improving safety.



Depositing box creates extra space for storage.



Foot switch is to operate the machine and dual-fuse safety design may shut off electricity when overloaded.



Standard Accessory

Suction tube (S)...180cm	2 pcs
Suction tube (L)...180cm	1 pc
Suction catheter	1 pc
Foot switch	1 pc
Filter (L)	2 pcs
Adaptor (D)	1 pc

Specification

Vacuum range	> 650 mmHg
Max. flow rate	120 L/min
Capacity	3000 c.c. + 3000 c.c.
Safety device	Floating control
N.W. / G.W.	24 kg / 28 kg
Machine size	L45 x W40 x H90 (cm)
Packing size	L52 x W48 x H94 (cm)

90-056 Oil-less Piston Vacuum Pump



Thermal protector, auto reset.

Pump specification

Voltage	120V / 240V
Current	3.5A / 1.7A
Capacitor	40uf / 13uf
R.P.M.	1720 / 1450 rpm
Output	375 W
Max. flow rate	120 L/min
Vacuum range	> 650 mmHg
Weight	8.3 kg



DOCTOR'S FRIEND medical instrument co., Ltd.

No.28, Ln. 13, Sec. 3, Zhongshan Rd., Tanzi Dist., Taichung City 42756, Taiwan
TEL : +886-425347777 FAX : +886-425338669 Email : medical.df@msa.hinet.net

